


## ANALYTICAL REPORT

## MMB-022 (MMB-4en-PICA) (C20H26N2O3)

## methyl 3-methyl-2-[[1-(pent-4-en-1-yl)-1H-indol-3-yl]formamido]butanoate

Remark – other NPS detected: **none**

Sample ID:	1955-18
Sample description:	crystallinic
Sample type:	test purchase /ISF projekt (NFL-SI)
Date of sample receipt (DD/MM/YYYY):	26/07/2018
Date of entry (DD/MM/YYYY) into NFL database:	13/08/2018
Report updates (if any) will be published here:	<a href="http://www.policija.si/apps/nfl_response_web/seznam.php">http://www.policija.si/apps/nfl_response_web/seznam.php</a>

Substance identified - structure <sup>1</sup> (base form)	
Systematic name	methyl 3-methyl-2-[[1-(pent-4-en-1-yl)-1H-indol-3-yl]formamido]butanoate
Other names	methyl 3-methyl-2-[[1-(pent-4-en-1-yl)-1H-indole-3-carbonyl]amino]butanoate; Methyl 3-methyl-2-(1-(pent-4-en-1-yl)-1H-indole-3-carboxamido)butanoate; methyl (1-(pent-4-en-1-yl)-1H-indole-3-carbonyl)valinate, MMB-022; MMB-4en-PICA;
Formula (per base form)	C20H26N2O3
M <sub>w</sub> (g/mol)	342,44
Salt form/anions detected	base
StdInChIKey (per base form)	RQVMOMGJHFPBNR-UHFFFAOYSA-N
Other NPS detected	none
Additional info (purity..)	> 96 % based on H1 NMR

<sup>1</sup> Created by OPSIN free tool: <http://opsin.ch.cam.ac.uk/> DOI: 10.1021/ci100384d

## Report updates

date	comments (explanation)

### Instrumental methods (if applied) in NFL

**1. GC-MS** (Agilent): GC-method is RT locked to tetracosane (9.258 min). Injection volume 1 ml and split mode (1:50). Injector temperature: 280 °C. Chromatographic separation: on column HP1-MS (100% dimethylpolysiloxane), length 30 m, internal diameter 0.25 mm, film thickness 0.25 µm. Carrier gas He: flow-rate 1.2 ml/min. GC oven program: 170 °C for 1 min, followed by heating up to 190 °C at rate 8 °C/min, then heating up to 293 °C at a rate of 18 °C/min, hold for 7.1 min, then heating at 50 °C/min up to 325 °C and finally 6.1 min isothermal. MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadropole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (30 until 6 min.) to 550 (300 until 6 min) amu.

**2. HPLC-TOF** (Agilent): 6230B TOF with Agilent 1260 Infinity HPLC with binary pump, column: Zorbax Eclipse XDB-C18, 50 x 4.6 mm, 1.8 micron. Mobile phases (A) 0.1% formic acid and 1mM ammonium formate in water; (B) 0.1% formic acid in methanol (B). Gradient: starting at 5% B, changing to 40% B over 4 min, then to 70% over 2 min and in 5 min to 100%, hold 1 min and back to 5%, equilibration for 1.7 min. The flow rate: 1.0 ml/min; Injection volume 1 µl. MS parameters: 2GHz, Extended Dynamic range mode to a maximum of 1700 amu, acquisition rate 1.30 spectra/sec. Sample ionisation: by Agilent Jet Stream technology (Dual AJS ESI). Ion source: positive ion scan mode with mass scanning from 82 to 1000 amu. Other TOF parameters: drying gas (N<sub>2</sub>) and sheath temperature 325 °C; drying gas flow rate 6 l/min; sheath gas flow rate 8 l/min; nebulizer 25 psig; Vcap. 4000 V; nozzle 2000 V; skimmer 65 V; fragmentor 175 V and Octopole RF 750 V.

**3. FTIR-ATR** (Perkin Elmer): scan range 4000-400 cm<sup>-1</sup>; resolution 4cm<sup>-1</sup>

**4. GC- (MS)-IR** condensed phase (GC-MS (Agilent) & IR (Spectra analyses-Danny)

GC-method: Injection volume 1 ml and split mode (1:5). Injector temperature 280 °C. Chromatographic separation as above (**1**). Split MS : IR = 1: 9.

MSD source EI = 70 eV. GC-MS transfer line T= 235°C, source and quadropole temperatures 280°C and 180°C, respectively. Scan range m/z scan range: from 50 (30 until 6 min.) to 550 (300) amu.

IR (condensed (solid) phase): IR scan range 4000 to 650, resolution 4 cm<sup>-1</sup>.

**5. IC** (anions) (Thermo Scientific, Dionex ICS 2100), Column: IonPac AS19, 2 x 250mm; Eluent: 10mM KOH from 0 to 10 min, 10-58 mM from 10 to 40min; Flow rate: 0.25 ml/min; Temperature: 30°C; Suppressor: AERS 500 2mm, suppressor current 13mA; Inj. Volume: 25 µl


## Supporting information

Solubility in	result/remark
CH <sub>2</sub> Cl <sub>2</sub>	soluble
MeOH	soluble
H <sub>2</sub> O	partially



Analytical technique:	applied	remarks
GC-MS (EI ionization)	+	NFL GC-RT (min): 11,1 BP(1): 212; BP(2): 227,BP(3) :144,
HPLC-TOF	+	Exact mass (theoretical): 342,1943; measured value Δppm:0,28; formula:C <sub>20</sub> H <sub>26</sub> N <sub>2</sub> O <sub>3</sub>
FTIR-ATR	+	direct measurement (sample as received)
FTIR (solid phase) always as base form	+	
IC (anions)	+	spot tests only
NMR (in FKKT)	+	
validation		
other		

# ANALYTICAL RESULTS



MS (EI 70 eV)



FTIR-ATR - direct measurement (sample as received)




IR (solid phase – after chromatographic separation)




University  
of Ljubljana  
Faculty of Chemistry  
and Chemical Technology



## R E P O R T

Contract No.	C1714-17-460078 (Republic of Slovenia, Ministry of the Interior, POLICE)
Sample ID:	<b>1955-18</b>
Received date:	August 6, 2018
Our notebook code:	NFL-1955-18
NMR sample preparation:	26.0 mg dissolved in 0.7 mL DMSO- <i>d</i> <sub>6</sub>
NMR experiments:	<sup>1</sup> H, <sup>13</sup> C, <sup>1</sup> H- <sup>1</sup> H <i>gs</i> -COSY, <sup>1</sup> H- <sup>13</sup> C <i>gs</i> -HSQC, <sup>1</sup> H- <sup>13</sup> C <i>gs</i> -HMBC, <sup>1</sup> H- <sup>15</sup> N <i>gs</i> -HMBC
Proposed structure with formula, exact mass, molecular weight:	 <p>Chemical Formula: C<sub>20</sub>H<sub>26</sub>N<sub>2</sub>O<sub>3</sub> Exact Mass: 342,1943 Molecular Weight: 342,4390</p>
Chemical name:	methyl (1-(pent-4-en-1-yl)-1 <i>H</i> -indole-3-carbonyl)valinate
Comments:	- Structure elucidation based on 1D and 2D NMR spectra and HRMS. - >96% purity of a sample based on <sup>1</sup> H NMR spectrum
Supporting information:	Copies of <sup>1</sup> H and <sup>13</sup> C NMR spectra, <sup>1</sup> H and <sup>13</sup> C FIDs.
Principal investigator:	Prof. Dr. Janez Košmrlj
Date of report:	August 10, 2018

NFL-1955-18  
1H



```


Current Data Parameters
NAME      NFL-1955-18
EXPNO    1
PROCNO   1

F2 - Acquisition Parameters
Date_    20180806
Time     23:03
INSTRUM  spect
PROBHD   5 mm PABBO BB-
PULPROG  zg30
TD       65536
SOLVENT  DMSO
NS       32
DS       2
SHE      10000.000 Hz
FIDRES   0.152588 Hz
AQ       3.2767999 sec
RG       57
DW       50.000 usec
DE       6.50 usec
TE       296.0 K
D1       1.00000000 sec
TD0      1

===== CHANNEL f1 =====
SFO1    500.1330885 MHz
NUC1    1H
P1      8.70 usec
PLM1    26.00000000 W

F2 - Processing parameters
SI      65536
SF      500.1300046 MHz
WDW     EM
SSB     0
GB      0
PC      1.00
  
```





Current Data Parameters  
 NAME NFL-1955-18  
 EXPNO 2  
 PROCNO 1

F2 - Acquisition Parameters  
 Date\_ 20180807  
 Time 1.47  
 INSTRUM spect  
 PROBHD 5 mm PABBO BB-  
 PULPROG zgpg30  
 TD 65536  
 SOLVENT DMSO  
 NS 3072  
 DS 4  
 SWE 29761.904 Hz  
 FIDRES 0.454131 Hz  
 AQ 1.1010048 sec  
 RG 2050  
 DW 16.800 usec  
 DE 6.50 usec  
 TE 296.0 K  
 D1 2.0000000 sec  
 D11 0.0300000 sec  
 TD0 1

==== CHANNEL F1 =====  
 SFO1 125.7703637 MHz  
 NUCL1 13C  
 P1 8.70 usec  
 PLW1 122.0000000 W

==== CHANNEL F2 =====  
 SFO2 500.1320005 MHz  
 NUCL2 1H  
 CPDPRG12 waltz16  
 PCPD2 80.00 usec  
 P1M2 26.00000000 W  
 P1M12 0.30046001 W  
 PLW13 0.15113001 W

F2 - Processing Parameters  
 SI 32768  
 SF 125.7578452 MHz  
 WDW EM  
 SSB 0  
 IB 1.00 Hz  
 GB 0  
 PC 1.40